





Kegel8 Tight & Tone Quick Start Guide



1. Insert the battery

Remove the back cover and insert 1 x 9 volt battery before replacing the cover. Don't forget to remove the black battery cap.



5. Insert your probe Squeeze a small amount of water-based lubricant onto the metal parts of the probe and insert vaginally. Ensure the electrodes point hip-to-hip.



2. Clean your probe

Give your probe a quick clean with a body-friendly, antibacterial spray such as Kegel8® Cleaning Spray.

ON

OF


6. Select a programme

Choose your programme by

pressing the PRG button. (See


the back page for help selecting

the right programme for you.)




3. Connect your probe

Insert the lead wire into the top of unit (either channel A or B) and connect your probe to the lead wire connectors. **Note:** some probes connect straight into the device and do not need a lead wire.



4. Turn on your Kegel8® Tight & Tone Press the on/off button once to turn on your device.



7. Adjust the power If using Channel A press the A+ button until you feel a slight tingle. You can decrease this using the A- button.



8. Exercise Start your programme and

begin enjoying the benefits of a stronger pelvic floor!

9. Clean and store

Clean your probe and store with your unit in the carry case provided.

Choosing the right programme for you

The Kegel8 Tight & Tone has a selection of programmes tailored to help you improve your pelvic floor strength.

Programme		Duration
P01	Pain Relief	20 minutes
P02	Urge (Sudden need to urinate)	20 minutes
P03	Stress Incontinence 1 (Mild)	20 minutes
P04	Stress Incontinence 2 (Severe)	20 minutes
P05	Frequency/Urge1	20 minutes
P06	Frequency/Urge 2 (Frequent & sudden urge to urinate even when the bladder is empty)	15 minutes
P07	Lack of Sensitivity	25 minutes
P08	Pelvic Floor Work Out	60 minutes
P09	Building Up Endurance	20 minutes

For the first week of your routine, increase the mA until you feel a slight tingling sensation. After this, increase your mA's 1-2 levels every 1-2 weeks. The toner can be used once a day.

Contact Us

Call us: +44 (0) 1482 496 932

Email us: info@kegel8.co.uk

Visit us: www.Kegel8.co.uk

Our friendly advisors are available Monday to Friday, 8am – 4pm (Closed bank holidays)

You can feel comfortable in asking us anything. We are not just a sales team. We are here to help you.