# Kegel 8




## Kegel8<sup>®</sup> Ultra 20 Quick Start Guide




#### 1. Insert the battery

Remove the back cover and insert one 9 volt battery before replacing the cover. Don't forget to remove the black battery cap.



#### 3. Connect your probe

Insert the lead wire into Channel A at the top left of your unit and connect your probe to the lead wire connectors **Note:** some probes connect straight into the device and do not need a lead wire.




#### 5. Insert your probe

Squeeze a small amount of water-based lubricant, such as Kegel8 Conductivitiy Gel on the metal parts of the probe and insert.



#### 2. Clean your probe

Give your probe a quick clean with a body-friendly, antibacterial spray such as Kegel8 Antibacterial Cleaning Spray.



#### 4. Turn on your Kegel8® Ultra 20

Press the on/off button once to turn on your device.



#### 6. Select a programme

Choose your programme by pressing the PRG button. See the seperate programme guide for more information.



#### 7. Exercise

Press the A+ button until you feel until you feel a slight tingle. You can decrease this using the A- button. Start your programme and enjoy the benefits of a stronger pelvic floor!

#### 8. Clean and store

Clean your probe and store in the provided carry case with your unit. It's best to store your probe with the lead wires still attached.

#### First timer? Let's test your sensation levels

If this is your first time using the Kegel8<sup>®</sup> Ultra 20 with a probe, then start with the **P01 Sensation test**. Increase the mAs until you feel a warm tingling sensation (No higher than 50mA).

If you feel a sensation, move onto PO2 Sensation Test.

If you feel **no** tingling sensation, perform **P01** at 50mA (or less) once a day for 4 weeks before moving onto **P02**.


## NOTE: If you cannot feel sensation returning after the 4 weeks, please consult a medical professional for advice.

**PO2 Sensation Test** - If you feel a tingling sensaiton below 50mA, progress to a treatment programme. If you cannot feel a sensation, use **PO2** once a day for a further 4 weeks.

#### Where to place your electrode pads

You can use electrode pads on the skin for a variety of Kegel8<sup>®</sup> Ultra 20 programmes. Simply connect your lead wire to the electrode pads and insert into channel A or B at the top of the unit. Use the corresponding controls to change the intensity. See below for examples of where you can place the pads and for more advice visit www.kegel8.co.uk







Abdomen

Sacrum

Ankle

### Choosing the right programme for you

The Kegel8<sup>®</sup> Ultra 20 has a selection of programmes tailored to help you improve your pelvic floor strength.



Muscle exercises - for strengthening your pelvic floor muscles

**Neuromodulation programmes -** for calming overactive nerves, circulation, relaxation and pain relief

Programme		Duration	Probe	Electrode Pads (recommended position)
P01	Sensation test 1	30 mins	~	✓ (your choice)
P02	Sensation test 2	60 mins	~	✓ (your choice)
P03	Pelvic Floor Muscle - Exercises 1	45 mins	~	✓ (sacrum)
P04	Pelvic Floor Muscle - Exercises 2	30 mins	<b>~</b>	✔ (sacrum)
P05	Pelvic Floor Muscle - Exercises 3 (Maintenance programme)	45 mins	~	✔ (sacrum)
P06	Sensation Improvement	20 mins	<b>~</b>	✔ (sacrum)
P07	Pelvic Floor Muscle - Exercises 4 (Endurance)	60 mins	~	✔ (sacrum)
P08	Sensory modulation (TENS)	60 mins	~	✓ (sacrum)
P09	Overactive Bladder neuro-modulation	20 mins	~	✓ (sacrum or ankle)
P10	Stress Incontinence 1	40 mins	~	✓ (sacrum)
P11	Stress Incontinence 2	40 mins	~	✓ (sacrum)
P12	Mixed Incontinence	30 mins	~	✔ (sacrum)
P13	Functional Continence Training	9 mins	~	✓ (sacrum)
P14	Postnatal / Post-Surgical Care 1	35 mins	~	✔ (sacrum)
P15	Postnatal / Post-Surgical Care 2	30 mins	$\checkmark$	✓ (sacrum)
P16	Bowel (Faecal) Incontinence 1	10 mins	✓	✔ (sacrum)
P17	Bowel (Faecal) Incontinence 2	40 mins	<b>v</b>	✔ (sacrum)
P18	Chronic Pain Relief	00	~	✓ (abdomen/sacrum)
P19	Pelvic Floor Circulation	60 mins	✓	✓ (abdomen/sacrum)
P20	Pelvic Floor Relaxation	20 mins	<b>v</b>	✔ (sacrum)

#### Please read the full manual for contraindications and further advice.

Not sure which<br/>programme is the<br/>right one for you?Call us: +44 (0) 1482 496 932<br/>(Monday to Friday, 8am - 4pm)

Email us: info@kegel8.co.uk Visit us: www.Kegel8.co.uk