Periform®+


Intra-Vaginal Probe

The Periform® has been the probe of choice for pelvic floor muscle stimulation and biofeedback for many of the worlds leading healthcare professionals for over thirteen years. After extensive user consultation, the Periform®+ is refined to improve the exceptional performance even further, whilst retaining the unique, patented and award winning form.



Unique Indicator Wand

- Visual aid to teach correct pelvic floor exercises
- New push fit
- New extension pieces provide greater amplification of pelvic floor movements to improve visual biofeedback
- Demonstrates benefit of bracing before coughing etc.
- Provides vital information on the efficiency of muscle stimulation






About Periform®+

- Can be used with most muscle stimulators and biofeedback equipment (refer to manufacturers instructions before use)
- Single patient use
- Medical grade surgical stainless steel
- Unique 'rectangular' section resists lateral movement
- Contoured electrodes provide comfortable stimulation

New features:

- New smooth clean shape
- Longer integrated cables
- New push fit indicator wand with extension pieces



Product Code 091253723